

The EsySuit™ Range



**LESS DRAPING,
MORE SURGERY**



EsySuit™ Supine



EsySuit™ Knee



EsySuit™ Lateral



EsySuit™ Shoulder



EsySuit™ Table

EsySuit™ Range

THE 1-MINUTE-DRAPING-SOLUTION

User-friendly draping in less than 1 minute with the EsySuit™ Range.

The time needed to perform the surgery itself is usually fairly constant. The time needed for preparing the operating room and the patient however is subject to variability and therefore less predictable.

A significant amount of OR time can be saved by standardizing the draping procedure with the EsySuit™ Range. Not only is the draping time reduced to less than 1 minute but also the disinfection of the patient may take less time. The EsySuit™ Range delivers a positive contribution to the use of resources and the associated efficient functioning of the OR.



EASE OF USE

Due to its unique patented design, EsySuit™ drapes can be applied in less than 1 minute.

Even when the OR occupancy is subject to variability, the application of the EsySuit™ Range is quick and easy to master.



EFFICIENT STAFF ALLOCATION

Compared to traditional draping, the number of handlings with the EsySuit™ drapes can be reduced to a minimum. As a result, different processes run parallel which creates an optimal OR Pulse.



INCREASED ERGONOMICS

The EsySuit™ Range is designed in such a way that nurses are less exposed to a non-ergonomic position while applying the EsySuit™ drapes.



TIME SAVING

The EsySuit™ Range is a one drape solution in which incision window and foil* are already integrated. This allows a standardized and reproducible way of draping which is a crucial part of an efficient and cost-effective OR.



ENVIRONMENTALLY FRIENDLY

The EsySuit™ Supine drape covers 9,6m² and generates approximately between 50% and 66% less waste than traditional draping solutions.¹ The EsySuit™ Range makes a positive contribution to the OR waste management.



IMPROVED EFFICIENCY

The efficient and safe functioning of the OR is paramount for every hospital. The EsySuit™ Range makes a positive contribution to this goal.

Item	Description	Art. code	Unit	Packed per
EsySuit™ Supine	One minute draping set supine	ME-ED 1000	PC	5
EsySuit™ Lateral	One minute draping set lateral	ME-ED 1100	PC	5
EsySuit™ Table	One minute draping set table	ME-ED 1200	PC	4
EsySuit™ Knee	One minute draping set knee	ME-ED 2000	PC	5
EsySuit™ Shoulder	One minute draping set shoulder	ME-ED 3000	PC	5

*Not applicable on EsySuit™ Knee

¹ MedEnvision Internal Data: Comparison of Operating Room Waste for Multi-piece Draping System vs. the EsySuit™ Range One Minute Draping System.